

# INI-43

Catalog No: #S08726



Package Size: #S08726-1 1 mg #S08726-2 10 mg #S08726-3 100 mg #S08726-4 200 mg #S08726-5 25 mg #S08726-6 5 mg #S08726-7 50 mg  
Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	INI-43
Other Names	AP-1;inhibit;Inhibitor;INI43;INI-43;INI 43;Kpnβ1;Kpnβ1 cargoes;NFAT;NFY;NFκB
Calculated MW	385.46
Storage	Powder: -20°C for 3 years . In solvent: -80°C for 1 year . Shipping with blue ice.

## Application Details

ε ε• εO• :C22H23N7Smile:CN(C)CCN1c(N)c(-c2nc3ccccc3[nH]2)c2nc3ccccc3nc12DMLε• :MFCD06650809η1—ε1;99.72%ζIAθ\$£ε1;DMSO:3.85 mg/mL  
(10 mM)δΩ'ε• ·ι θ—:ε 'Aδ» ι• Aη 'H:Others

Note: This product is for in vitro research use only and is not intended for use in humans or animals.