

ISAM-140

Catalog No: #S08734



Package Size: #S08734-1 1 mg #S08734-2 10 mg #S08734-3 100 mg #S08734-4 200 mg #S08734-5 25 mg #S08734-6 5 mg #S08734-7 50 mg
Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	ISAM-140
Other Names	ISAM 140;ISAM140;ISAM-140;A2B;Adenosine Receptor;AdenosineReceptor
Calculated MW	337.37
Storage	Powder: -20°C for 3 years . In solvent: -80°C for 1 year . Shipping with blue ice.

Application Details

SMILES: C19H19N3O3
 Molecular Weight: 337.37 mg/mL (10 mM)
 A2B (human): 3.49 nM (Ki)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.