

# Panaxadiol

Catalog No: #S08735



Package Size: #S08735-1 1 mL #S08735-2 100 mg #S08735-3 200 mg #S08735-4 25 mg #S08735-5 50 mg  
 Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
 Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Panaxadiol
Other Names	anti-arrhythmic activity;antioxidative activity;anticancer activity;20(R)-Panaxadiol;20Panaxadiol;cardioprotective activity;HIF;HIF ProlylHydroxylase;HIF Prolyl-Hydroxylase;HIF/HIFProlylHydroxylase;HIF/HIF ProlylHydroxylase;HIFProlylHydroxylase;Panaxadiol;SW620;δ1ε• δ1 i
Calculated MW	460.73
Storage	keep away from direct sunlight,keep away from moisture . Powder: -20°C for 3 years . In solvent: -80°C for 1 year . Shipping with blue ice.

## Application Details

ε ε• εO• :C30H52O3Smile:[H][C@@]1(CC[C@]2(C)[C@]1([H])[C@H](O)C[C@]1([H])[C@@]3(C)CC[C@H](O)C(C)(C)[C@]3([H])CC[C@@]21C)[C@@]1(C)CCCC(C)(C)O1DMLε• :MFC00210540ηl—ε1;99.9%ζIAθ§£ε1;Ethanol:20 mg/mL (43.41 mM)δΩ'ε• ·i θ·—:θ"θ§ i• δO |δ»£θ"i• Aη 'H:HIF/HIF Prolyl-Hydroxylase

Note: This product is for in vitro research use only and is not intended for use in humans or animals.