

Isbogrel

Catalog No: #S08736



Package Size: #S08736-1 1 mg #S08736-2 10 mg #S08736-3 100 mg #S08736-4 25 mg #S08736-5 5 mg #S08736-6 50 mg
Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

| | |
|---------------|---|
| Product Name | Isbogrel |
| Other Names | CV 4151;CV-4151;CV4151;Isbogrel;Prostaglandin Receptor;ProstaglandinReceptor;δO ζ³ζ Oι ;Thromboxane (TX) synthase |
| Calculated MW | 281.35 |
| Storage | Powder: -20°C for 3 years . In solvent: -80°C for 1 year . Shipping with blue ice. |

Application Details

ε ε• εO• :C18H19NO2Smile:C(=C\C\CCCC(O)=O)\(C1=CC=CC=C1)/C=2C=CC=NC2DMLε• :MFCDD01753219η!—ε!;97.06%ζ!Aθ§εε!;DMSO:55 mg/mL (195.49 mM)δ½ ε ζ»ζ §:Isbogrelι θΩ ε• Εζ • (1ε 10 mg/kg)ε ι• θ ζ³ε° (1 mg/kg)ζ `Υθ ζ ε Αδ! ε ε ε-|θ—±ε—Oη MCAθ' ζ ε½ζ • γ [2] ζ- ε οO Isbogrel (10 mg/kg, p.o.)ζ `Υθ ε • ε°• δ! ζ' ε' ε 1ε οO ζ ε Αδ! δH³ι Έε• «ι • η ε' ε γ ε ζ- οO Isbogrelε• —θ ½ζ —ζ²»η θ θ' ζ ε ι' ι ²θ ζ' ε' η ζ ζ ε ε• η ©γ [2]δΩ'ε• ι θ—:G θ η ½ε• Αθ• ε• δ½ |ε • η «δE η η ι• Άη 'H:Prostaglandin Receptor

Note: This product is for in vitro research use only and is not intended for use in humans or animals.