

Ispronicline

Catalog No: #S08740



Package Size: #S08740-1 1 mL

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Ispronicline
Other Names	RJR-1734;TC 01734;TC 1734;TC-1734;TC1734;TC-01734;TC01734;Ispronicline;AZD3480;AZD-3480;AZD 3480;εO δE ε ε °;α4β2 nAChR
Calculated MW	234.34
Storage	Powder: -20°C for 3 years . In solvent: -80°C for 1 year . Shipping with blue ice.

Application Details

ε ε• εO• :C14H22N2OSmile:CN[C@@H](C)C\C=C\Cc1cncc(OC(C)C)c1DMLε• :MFCD11040931ηl—εl;96.54%ζlAθ§£εl;DMSO:50 mg/mL (213.37 mM)δ½ ε ζ»ζ §:ε °ε §lO δE-oO η δE ζ°ε° IsproniclineoO 10-20 mg/kgoO θ ½ε ε °ε • θ η η Hε ε lε εO ε• c-Fosθ°θY°Yη θ—±ε—O[1]γ Isproniclineη θ• —δ»£ε °ε ε-η HεY• ζ —ε °ε §lO δ½ ε ε• θ°ζ δE12ε° ζ AoO δE ζ- ε½ζ • ε—Hζ— η ζ —oO ε θ ½ε ζ • η»-ε' εO1ε• δ½ θ °εΩ oO ζ • η»-ζ Aι 'δ» 18ε° ζ Ae »Aη»-θ °2ε ©[2]γ δΩ'ε• ι θ—:η η»• η§ ε-ι• Aη H:AChR

Note: This product is for in vitro research use only and is not intended for use in humans or animals.