

# JAK3i

Catalog No: #S08745



Package Size: #S08745-1 1 mg #S08745-2 1 mL #S08745-3 10 mg #S08745-4 25 mg #S08745-5 5 mg #S08745-6 50 mg  
Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	JAK3i
Other Names	JAK 3i;JAK3;JAK3i;JAK-3i
Calculated MW	354.34
Storage	keep away from moisture . Powder: -20°C for 3 years . In solvent: -80°C for 1 year . Shipping with blue ice.

## Application Details

ε ε• εO• :C18H15FN4O3Smile:C(OCC)(=O)C=1C2=C(N=CN=C2NC1)C3=CC(NC(C=C)=O)=CC=C3Fη1—ε1!;99.81%ζ1Aθ\$εε1!;DMSO:10 mg/mL (28.22 mM)δ½ ε ζ'»ζ §:ζ δ»γ ε—Hε ρ η δE ζ'E'IL2Rδ» δ½ JAK

3οO δ½Ωη "δI δE η§• δE'εI η 'Eε ρη ι« ε!;η HεO ζ §JAK3ζ ε 'Aε (JAK3i; PF-06651600)οO θ— ζ ε Aε ε "δ½ ε εO1ζ ζ ε AδI STAT5η ηε·I 'Eε γ δ½ ε ι • JAK3iδ'E η» θ ε θ ½η «θ η- η η 'Eη» ε• οO η 'Eζ— δI ε• η ρη δ»»δ'E η ζ³ οO θΩ δ'E ζ- ι • δ½ δI θ Ωη θ' θ• γ [1]

ι • ε• η θ"ε ζ ρδ½ η» ζ ζ• -η 'δI JAK3ιγ • PLCG2ε ZEB2δ½ δ'E1ε ρι ζ 'η1½ε 1ε οO ε δ»-ε—HδI δ'E εΩ θ' η 'ε ε • η «ε 1 η- η 'Eε ρθ"ε η θ ε• η η ζ η• θ"η γ [2]δΩ'ε• ι θ—:θ' η 'η ζ • |θ"θ§ ι• δO |JAK/STAT δΩ'ε• ι θ—|θ η ½ι ζ"ι 'EζΩ ι A|εH²η» θ ι• Aη H:JAK

Note: This product is for in vitro research use only and is not intended for use in humans or animals.