

# JNJ-39220675

Catalog No: #S08754



Package Size: #S08754-1 1 mg #S08754-2 10 mg #S08754-3 25 mg #S08754-4 5 mg #S08754-5 50 mg  
Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	JNJ-39220675
Other Names	JNJ 39220675;JNJ39220675;JNJ-39220675;HistamineReceptor;H3 receptor;Histamine Receptor
Calculated MW	369.43
Storage	Powder: -20°C for 3 years . In solvent: -80°C for 1 year . Shipping with blue ice.

## Application Details

O=C(C1=CN=C(OC2=CC=C(F)C=C2)C=C1)N3CCN(CCC3)C4CCC4DML
  
 Molecular Weight: 369.43
   
 Purity: >99.78%
   
 Storage: Powder: -20°C for 3 years . In solvent: -80°C for 1 year . Shipping with blue ice.
   
 Target: Histamine Receptor

Note: This product is for in vitro research use only and is not intended for use in humans or animals.