

JNJ-42253432

Catalog No: #S08759



Package Size: #S08759-1 1 mg #S08759-2 1 mL #S08759-3 10 mg #S08759-4 100 mg #S08759-5 200 mg #S08759-6 2.5 mg #S08759-7 5 mg #S08759-8 50 mg
Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	JNJ-42253432
Other Names	JNJ 42253432; JNJ42253432; JNJ-42253432; P2X Receptor; P2X7; P2X Receptor
Calculated MW	446.63
Storage	Powder: -20°C for 3 years . In solvent: -80°C for 1 year . Shipping with blue ice.

Application Details

ε ε• εO• :C28H38N4O
 SMILES: CN1CCc2c(C1)ccc2C(=O)NCC1(CCCCC1)N1CCN(CC1)c1cccc1η1—ε1; 99.08% ζ|Aθ\$£ε1; DMSO: 90.0 mg/mL (201.5 mM) • Άη 'Hζ' »ζ §: P2X7: 9.1 (pKi, rat) | P2X7: 7.9 (pKi, human) δΩ'ε• ·ι θ·—: η|»ε• ι ι• |η η»• η§ ε-|ι• Άη 'H: P2X Receptor

Note: This product is for in vitro research use only and is not intended for use in humans or animals.