

3-Aminopropionitrile fumarate (2:1)

Catalog No: #S08760



Package Size: #S08760-1 1 g #S08760-2 1 mL #S08760-3 10 g #S08760-4 100 mg #S08760-5 5 g #S08760-6 500 mg
 Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	3-Aminopropionitrile fumarate (2:1)
Other Names	3-Aminopropionitrile fumarate;3Aminopropionitrile fumarate (2:1);3-Aminopropionitrile fumarate (2:1);3-Aminopropionitrile fumarate 2:1;3 Aminopropionitrile fumarate (2:1);2:1;beta-Ammoniumpropionitrile;beta-Ammoniumpropionitrile Hemifumarate;Beta-Aminopropionitrile fumarate;b-Aminopropionitrile Fumarate;b-Ammoniumpropionitrile;b-Ammoniumpropionitrile Hemifumarate;Di-β-aminopropionitrile fumarate;Di-beta-aminopropionitrile Fumarate;Di-b-aminopropionitrile Fumarate;ε— i@~ı 'E3-ζ''ε 1δ'E θ ;β-Aminopropionitrile;β-Aminopropionitrile fumarate;β-Ammoniumpropionitrile;β-Ammoniumpropionitrile hemifumarate
Calculated MW	114.12
Storage	Powder: -20°C for 3 years . In solvent: -80°C for 1 year . Shipping with blue ice.

Application Details

ε ε• εO• :C3H6NB·1/2C4H4O4Smile:NCCC#N.NCCC#N.OC(=O)\C=C\O)=ODMLε• :MFCD00013142η1—ε1';99%ζ\Aθ\$£ε1';DMSO:12 mg/mL (46.83 mM)δΩ'ε• ·ı θ·—:ε Aδ» ı• Aη 'H:Others

Note: This product is for in vitro research use only and is not intended for use in humans or animals.