

JNJ-DGAT2-A

Catalog No: #S08761



Package Size: #S08761-1 1 mg #S08761-2 10 mg #S08761-3 100 mg #S08761-4 25 mg #S08761-5 5 mg #S08761-6 50 mg #S08761-7 500 mg
Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	JNJ-DGAT2-A
Other Names	JNJDGAT2A;JNJ-DGAT-2-A;JNJ-DGAT2-A;JNJ DGAT2 A;DGAT2
Calculated MW	523.38
Storage	Powder: -20°C for 3 years . In solvent: -80°C for 1 year . Shipping with blue ice.

Application Details

ε ε• εO• :C24H16BrFN4O2SSmile:Fc1cc(\C=C2/SC(NCCc3ccccc3)=NC2=O)ccc1Oc1ccc(cc1Br)C#NDMLε• ·:MFCD30182366η!—ε!;:98.2%ζ!Aθ§εε!;DM
SO:3.0 mg/mL (5.7 mM)δ½ ε ζ»ζ §:JNJ-DGAT2-A (5 ΞOM)
ε δH ε ε ζ ε Aδ! ι • η» DGAT2ε HepG2η» θ θ£ θ§£η ©δ'E-η DGAT2ζ'»ζ §γ ε "HepG2η» θ δ'E-oO JNJ-DGAT2-A (0.3125-20
μM)ζ ε ι • δY• θ" η ζ HεO• ζ ε Aδ! TG (50:2)γ • TG (52:2)ε TG (54:3)η η ζ • oO ε AIC50ε Oε ε «δE!0.66γ • 0.85ε 0.99
ΞOM[1]γ ι• Aη 'Hζ'»ζ §:DGAT-2 (Sf9 insect cell membranes):0.14 ΞOMδΩ'ε• ι θ·—:δ»£θ°ι• Aη 'H:Transferase

Note: This product is for in vitro research use only and is not intended for use in humans or animals.