

JTK-109

Catalog No: #S08763



Package Size: #S08763-1 1 mg #S08763-2 10 mg #S08763-3 100 mg #S08763-4 25 mg #S08763-5 5 mg #S08763-6 50 mg
Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	JTK-109
Other Names	HCV Protease;HCVProtease;JTK 109;JTK109;JTK-109;NS5B polymerase
Calculated MW	638.13
Storage	Powder: -20°C for 3 years . In solvent: -80°C for 1 year . Shipping with blue ice.

Application Details

ε ε• εO• :C37H33ClFN3O4Smile:FC1=C(C=2N(C=3C(N2)=CC(C(O)=O)=CC3)C4CCCCC4)C=CC(OCC5=C(C=CC(=C5)N6C(=O)CCC6)C7=CC=C(CI)C=C7)=C1DMLε• :MFCD16294956η!—ε!;98.68%ζ!Aθ§εε!;DMSO:6.38 mg/mL (10 mM)δΩ'ε• ι θ·—:εY η η ©ε-||θ η ½ Aδ½ ι• Aη 'H:HCV Protease

Note: This product is for in vitro research use only and is not intended for use in humans or animals.