

K-111

Catalog No: #S08767



Package Size: #S08767-1 1 mg #S08767-2 10 mg #S08767-3 25 mg #S08767-4 5 mg #S08767-5 50 mg
Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	K-111
Other Names	BM 170744;BM170744;BM-170744;K 111;K111;K-111;PPARα
Calculated MW	379.75
Storage	Powder: -20°C for 3 years . In solvent: -80°C for 1 year . Shipping with blue ice.

Application Details

ε ε• εO• :C18H25Cl3O2Smile:O=C(O)C(Cl)(Cl)CCCCCCCCC1=CC=C(Cl)C=C1η1—ε1!;99.93%δ½ ε ζ'»ζ §:K-111ζ²»η οO ζ—• 3ε 4ε ``ζ 'ζ
 §ε ι • οO 1οO 3ε 10mg/kg/dοO ε• οO K-111ε—Oθ 'δ½ ι • ζ 'Υθ δE ι • οO ζ 'Hε ι« θ °ε² η'θ' η γ • θ °ε² η'ζ • ζ ζ §Υ •
 ι« η ζ²HδE ι —θ' η ε ι« ε— ε1!θ θ η ½θ ε 'ιι ζ°εH°οO ε—'Hη©'ιθ 'Hθ' η³ γ • 24ε°• ζ Ae°QζA²ζ ζ³ γ • ζ AηO©ε ζ θ
 εOε γ • θ' ζ" η1 η»θ η ½ε γ • ζ »θ ε 'ιι ζ ε ε-|ε θ' ζA²ε-|ζ| ε "ζ²ζ ζ 'Υθ ε½±ε • γ [6]δΩ'ε• ·ι θ·—:δ»£θ°|DNA
 ζ• δO ε δΩ ε • ι• Aη 'H:PPAR

Note: This product is for in vitro research use only and is not intended for use in humans or animals.