

# K67

Catalog No: #S08768

Package Size: #S08768-1 2 mg



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	K67
Other Names	K67;K-67;K 67;Keap1-Nrf2
Calculated MW	582.69
Storage	store at low temperature . Powder: -20°C for 3 years . In solvent: -80°C for 1 year . Shipping with blue ice.

## Application Details

ε ε• εO• :C29H30N2O7S2Smile:O=C(C)CC=1C=C(NS(=O)(=O)C2=CC=C(OCC)C=C2)C=3C=CC=CC3C1NS(=O)(=O)C4=CC=C(OCC)C=C4DMLε• :MF CD31697709η1—ε1;98.43%ζ1Aθ§ε1;DMSO:80 mg/mL (137.29 mM)δ½ ε ζ»ζ §:ζ 'Hζ³ οO K67οO 50

ΞOMοO 4ε° ζ AοO ε η• θ"θ'YYS351Eη Huh7η» θ οO η» θ η ζ δ'E• η "δH³θ Ωζ± οO 10ΞOMοO ζ²»η 8ε° ζ Ay ε—HεEζ ζ Nrf2ζ δ½ η ε •η «ζ ε Ae θΩ θ' δ1 ε •η «ε• °θΩHe ζ •γ

η» ζ οO η "K67ε η• Huh1η» θ ζ •ε §ε °ζ ε Aδ1 p62η£·ι Έε ε Keap1δ'H ι 'η η Έδ1 δ½ η "οO K67ζ 'Yθ δΩ θΩ δ1 Nr f2η ζ³ η'ε γ [1]δΩ'ε• ·ι θ—:ε •η «δ'E η η ι• Aη 'H:Nrf2

Note: This product is for in vitro research use only and is not intended for use in humans or animals.