

K6PC-5

Catalog No: #S08770



Package Size: #S08770-1 1 mg #S08770-2 1 mL #S08770-3 10 mg #S08770-4 100 mg #S08770-5 200 mg #S08770-6 2.5 mg #S08770-7 5 mg #S08770-8 50 mg

Support: tech@signalwayantibody.com

Description

Product Name	K6PC-5
Other Names	K6PC 5;K6PC5;K6PC-5;SphK1;S1PReceptor;S1P Receptor
Calculated MW	343.5
Storage	Powder: -20°C for 3 years . In solvent: -80°C for 1 year . Shipping with blue ice.

Application Details

ε ε• εO• :C19H37NO4Smile:CCCCCCCC(=O)C(CCCCC)C(=O)NC(CO)CODMLε• :MFCD30480923ηl—εl;98.55%ζlAθ§εεl;DMSO:90 mg/mL (262 mM)δ½ ε ζ'»ζ §:ε ζ -θ'''δε θ • ε η ζ ζ— ε° iO oO 56ε "ε²• oO δε-oO 1%η K6PC-5θ'«ε± ι "εl η "2ε "γ θΩ η§• K6PC-5ε η• ζ 'Yθ ε' ε δl η η ζ •ηl η»'η» θ η ζ 'ι • ε θ Ae θ η ½η δl§η γ ε ζ- oO η η ε εl;δH ζ 'Yθ ε' ε [2]γ δ½ ε ζ'»ζ §:K6PC-5oO 1-10

ΞOMoO 24ε°• ζ AoO ε "ε ι • δY• θ" ζ §ζ 'HeO• δε oO ζ• ι« δl ζ-εε'Eεδl'θ""η θ§ θ""ε½ζ • η» θ oO NHEKsoO δε-η θ§ θ η ½ε η'§ε— θ η ½ζ°'εH³γ K6PC-5ι θΩ η» θ ε Ca2+δΩ'ε• ·δO ε—OδΩ θΩ θ§ θ""ε½ζ • η» θ η ε ε ε ε' ζ γ ζ- ε ο O K6PC-5ζΩ ζ'»p42/44η» θ ε δΩ'ε• ·θ° θ ζΩ ι Ae c-Junζ «η«—ζΩ ι Aη ηε·ι 'Eε [1]γ K6PC-5oO 1-10

ΞOMoO 24ε°• ζ AoO δΩ θΩ δl'ζ • ηl η»'η» θ η ε' ζ ε θ Ae θ η ½η ε• ζ • γ K6PC-5ε "δl'ζ • ηl η»'η» θ δε-εO θ"·η » θ ε Ca2+ζ" εl;oO [Ca2+]ioO η ζ —θ• '[2]γ K6PC-5oO 10, 25ε 50

ΞOMoO 48ε°• ζ AoO ζ 'Yη• ε • θ'½»EBOVζ ζ EA.hy926η» θ ζ θ 'η ζ ζ γ K6PC-5ζ 'Yθ ι • δ½ δl ζ ζ η» θ δε ζE ζA°δε-η η ζ— ζ»'εl;oO ε'HAδ» ε ι • δY• θ" η ζ 'HeO• ζ 'Yθ ε • ε° δl VP40η ε• «ι • [3]γ δΩ'ε• ι θ·—:G θ η ½ε• Aθ• ε• δ½ ι• Aη 'H:S1P Receptor

Note: This product is for in vitro research use only and is not intended for use in humans or animals.