

# KD-026

Catalog No: #S08771



Package Size: #S08771-1 1 mg #S08771-2 1 mL #S08771-3 10 mg #S08771-4 100 mg #S08771-5 200 mg #S08771-6 2.5 mg #S08771-7 5 mg #S08771-8 50 mg

Order: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
 Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	KD-026
Other Names	microsomal triglyceride transfer protein;KD-026;KD026;KD 026;Inhibitor;inhibit;MTP inhibitor;MTP;enterocytic-specific;dyslipidemia;SLx 4090;SLx4090;SLx-4090
Calculated MW	546.54
Storage	Powder: -20°C for 3 years . In solvent: -80°C for 1 year . Shipping with blue ice.

## Application Details

ε ε• εO• :C31H25F3N2O4S  
 smile:COc1cccc(C(=O)Nc2ccc3CN(CCc3c2)C(=O)Oc2ccccc2)c1-c1ccc(cc1)C(F)(F)FDMLE• :MFCD28502203η1—ε1!;99.73%ζ1  
 Aθ\$εε1!;DMSO:55 mg/mL (100.63 mM)δ½ ε ζ'»ζ §:ε ε §iO δ'E-oO KD-026oO 30 mg/kgοO ι θΩ ι • δ½ 50%η ι • ε• θ θ οO ED50 =  
 7 mg/kgοO [1]γ δ½ ε ζ'»ζ §:KD-026 ε—Hθ½½θ θ η ½Bζ 'Υη 1ε 1ζ ε Aδ½ η οO IC50 = 9.6  
 nMoO οO δ½ ε—HCaco-2η» θ δE-θ½½θ θ η ½A1η ε ζ³ ζ ε½±ε • [1]γ ι• Aη 'Hζ'»ζ §:MTP:8.0 nM|ApoB:9.6  
 nMδΩ'ε• ι θ—:δ»εθ°ι• Aη 'H:CETP

Note: This product is for in vitro research use only and is not intended for use in humans or animals.