

# KD-027

Catalog No: #S08773



Package Size: #S08773-1 1 mg #S08773-2 10 mg #S08773-3 100 mg #S08773-4 25 mg #S08773-5 5 mg #S08773-6 50 mg  
Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|               |  |
|---------------|--|
| Product Name  | KD-027   |
| Other Names   | KD 027;KD027;KD-027;SLx 2101;SLx2101;SLx-2101                                      |
| Calculated MW | 561.65   |
| Storage       | Powder: -20°C for 3 years . In solvent: -80°C for 1 year . Shipping with blue ice. |

## Application Details

ε ε• εO• :C25H35N7O6SSmile:O=C1N=C(NN2C1=C(N=C2CCC)C)C=3C=C(C=CC3OCC)S(=O)(=O)N(C)CC(=O)N(C)CC(=O)N(C)Cη1—ε1!;99.55%δΩ'ε• ·  
I θ—:δ»£θ°!• Aη 'H:PDE

Note: This product is for in vitro research use only and is not intended for use in humans or animals.