

KDM5-C49 HCl

Catalog No: #S08774



Package Size: #S08774-1 1 mg #S08774-2 1 mL #S08774-3 10 mg #S08774-4 100 mg #S08774-5 25 mg #S08774-6 5 mg #S08774-7 50 mg #S08774-8 500 mg
 Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	KDM5-C49 HCl
Other Names	HistoneDemethylase;Histone Demethylase;KDOAM 20 Hydrochloride;KDOAM20 Hydrochloride;KDOAM-20 Hydrochloride;KDOAM-20 hydrochloride;KDM5-C49 HCl(1596348-16-1 Free base);KDM5-C49 HCl;KDM5-C49 HCl(1596348-16-1 Free base);KDM5
Calculated MW	344.84
Storage	Pure form: -20°C for 3 years . In solvent: -80°C for 1 year . Shipping with blue ice.

Application Details

SMILES: Cl.CCN(CCN(C)C)C(=O)CNCc1cc(ccn1)C(O)=O
 DMSO: 9 mg/mL (26.1 mM)
 KDM5B: 160 nM | KDM5A: 40 nM | KDM5C: 100 nM
 Histone Demethylase

Note: This product is for in vitro research use only and is not intended for use in humans or animals.