

KGP94

Catalog No: #S08776



Package Size: #S08776-1 1 mg #S08776-2 1 mL #S08776-3 10 mg #S08776-4 100 mg #S08776-5 200 mg #S08776-6 25 mg #S08776-7 5 mg #S08776-8 50 mg

Order: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

| | |
|---------------|---|
| Product Name | KGP94 |
| Other Names | Inhibitor;inhibit;M2 marker;KGP 94;KGP94;KGP-94;Raw264.7 macrophage;prostate (PC-3ML);Metastasis;metastatic carcinoma;Cysteine Protease;CysteineProtease;Cathepsin;breast cancer cell;Cathepsin L |
| Calculated MW | 350.23 |
| Storage | Powder: -20°C for 3 years . In solvent: -80°C for 1 year . Shipping with blue ice. |

Application Details

ε ε• εO• :C14H12BrN3OSSmile:NC(=S)N\N=C(lc1cccc(O)c1)c1cccc(Br)c1ηl—εl;98.52%ζlAθ§£εl;DMSO:60 mg/mL (171.32 mM)δ½ ε ζ'»ζ §:ζ °ε η ε• 'η 'Hζ ι AL/Kζ ε Aε KGP94θ ½ε ε "δ½ ε ζ ε AM2ε ε·"ε ρη» θ η δY"θ'-oO ε'HAι • δ½ ε· "ε ρη» θ ε lζΩ δE η 4T1ε°• iO δH°θ 'η η» θ δY"θ'-γ KGP94ε η• θΩ ι • δ½ δl ε η§• δE M2η 'Eε ρη ζ εΩ η ©η θ"θY 'YoO θ"ζ ε• 'η 'Hζ ι ALη ζ'»ζ §ε• —θ ½ε—HIL-4i©±ε "η M0ε °M2η ε ε θ °ε ρι • θl• [1]γ δΩ'ε• ι θ—:θ η ½l Aδ½ ι• Aη 'H:Cysteine Protease

Note: This product is for in vitro research use only and is not intended for use in humans or animals.