

L-162,313

Catalog No: #S08782



Package Size: #S08782-1 1 mg #S08782-2 10 mg #S08782-3 100 mg #S08782-4 25 mg #S08782-5 5 mg #S08782-6 50 mg #S08782-7 500 mg
Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	L-162,313
Other Names	L162,313;L-162,313;L162313;L-162313;L 162,313;L 162313;angiotensin II receptor
Calculated MW	582.78
Storage	Powder: -20°C for 3 years . In solvent: -80°C for 1 year . Shipping with blue ice.

Application Details

SMILES: C30H38N4O4S2
 DMSO: 50 mg/mL (85.8 mM)
 RAAS

Note: This product is for in vitro research use only and is not intended for use in humans or animals.