

L-750,667 Trihydrochloride

Catalog No: #S08785



Package Size: #S08785-1 1 mg #S08785-2 10 mg #S08785-3 100 mg #S08785-4 200 mg #S08785-5 25 mg #S08785-6 5 mg #S08785-7 50 mg
 Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	L-750,667 Trihydrochloride
Other Names	L750,667;L-750667 triHydrochloride;L-750,667;L750,667 Trihydrochloride;L-750,667 triHydrochloride;L750,667 TriHydrochloride;L 750,667
Calculated MW	527.66
Storage	Powder: -20°C for 3 years . In solvent: -80°C for 1 year . Shipping with blue ice.

Application Details

SMILES: CC1=CC=C(N2CCN(CC2)CC3=CNC4=C3C=CC=N4)C=C1.Cl.Cl.Cl
 Molecular Weight: 527.66
 DMSO: 5 mg/mL (9.48 mM)
 D4 Receptor: 0.51 nM (Ki) | D4 Receptor: 0.16 nM (Kd)
 Dopamine Receptor

Note: This product is for in vitro research use only and is not intended for use in humans or animals.