

Lapisteride

Catalog No: #S08793



Package Size: #S08793-1 1 mg #S08793-2 10 mg #S08793-3 100 mg #S08793-4 25 mg #S08793-5 5 mg #S08793-6 50 mg #S08793-7 500 mg
Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Lapisteride
Other Names	CS 891;CS891B;CS-891;CS891;CS-891B;Lapisteride;SRD5A1;ζ ε Ηι θ 1
Calculated MW	464.64
Storage	Powder: -20°C for 3 years . In solvent: -80°C for 1 year . Shipping with blue ice.

Application Details

ε ε• εO• :C29H40N2O3Smile:C[C@@]12[C@]([C@]3([C@@]([C@]4(C)[C@@](CC3)(NC(=O)C=C4)[H])(CC1)[H])[H])(CC[C@@H]2C(NC(C)(C)C5=CC=C(OC)C=C5)=O)[H]DMLε• :MFCD09837691η1—ε1!;98.77%ζ1Aθ§£ε1!;DMSO:50 mg/mL (107.61 mM)δΩ'ε• ·ι θ·—:ε 'Aδ» ι• Aη 'H:Others

Note: This product is for in vitro research use only and is not intended for use in humans or animals.