

# LAS38096

Catalog No: #S08795



Package Size: #S08795-1 1 mg #S08795-2 10 mg #S08795-3 25 mg #S08795-4 5 mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	LAS38096
Other Names	LAS 38096;LAS-38096;LAS38096;AdenosineReceptor;Adenosine Receptor;A2B
Calculated MW	316.32
Storage	Powder: -20°C for 3 years . In solvent: -80°C for 1 year . Shipping with blue ice.

## Application Details

ε ε• εO• :C17H12N6O5  
 smile:N(c1cccnc1)c1ncc(-c2ccncn2)c(n1)-c1cccc1DMLε• :MFCD19703267  
 η!—ε!;98.71%ζ!Aθ\$εε!;DMSO:55 mg/mL (173.87 mM)|• Aη 'Hζ'»ζ §:A2A:>2500 nM (Ki)|A1:>1000 nM (Ki)|A3:>1000 nM (Ki)|A2B:17 nM (Ki)  
 δΩ'ε• ·ι θ•—:G  
 θ η ½ε• Aθ• ε• δ½ |η η»• η§ ε-|• Aη 'H:Adenosine Receptor

Note: This product is for in vitro research use only and is not intended for use in humans or animals.