

# Sanguinarine

Catalog No: #S08798



Package Size: #S08798-1 1 mg #S08798-2 10 mg #S08798-3 100 mg #S08798-4 25 mg #S08798-5 5 mg #S08798-6 50 mg  
 Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
 Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Sanguinarine
Other Names	Antioxidant;Apoptosis;Autophagy;Inhibitor;inhibit;Sanguinarium;Sanguinarin;Sanguinarine;Pseudochelerythrine;θ' ζ'Hη'±
Calculated MW	332.33
Storage	Powder: -20°C for 3 years . In solvent: -80°C for 1 year . Shipping with blue ice.

## Application Details

ε ε• εO• :C20H14NO4Smile:C[N+]=1C2=C(C=3C(C1)=C4C(=CC3)OCO4)C=CC=5C2=CC6=C(C5)OCO6DMLε• :MFCD00064925η|—ε|':99.18%ζ'1Aθ§εε  
 '!:DMSO:3.32 mg/mL (9.99 mM)δQ'ε• ·ι θ—:ζ°§ε θΩ ε |θ ε -|ε δ'ι• Aη 'H:Apoptosis|Antioxidant|Autophagy

Note: This product is for in vitro research use only and is not intended for use in humans or animals.