

Lirequinil

Catalog No: #S08810



Package Size: #S08810-1 1 mg #S08810-2 10 mg #S08810-3 100 mg #S08810-4 25 mg #S08810-5 5 mg #S08810-6 50 mg #S08810-7 500 mg
 Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	Lirequinil
Other Names	RO 413696;Ro 41-3696;RO413696;RO-413696;Ro-41-3696;RO41-3696;Lirequinil;GABAA;GABARceptor;GABA Receptor;ε @η ε Ηε°O
Calculated MW	448.94
Storage	Powder: -20°C for 3 years . In solvent: -80°C for 1 year . Shipping with blue ice.

Application Details

ε ε• εO• :C26H25CIN2O3Smile:C(=O)(C1=C2C=3C(CCN2C(=O)C(=C1)C4=CC=CC=C4)=CC=C(CI)C3)N5C[C@@H](OCC)CC5DMLε• :MFCD00868037
 ηI—εI!;98.06%ζIAθ\$εI!;DMSO:50 mg/mL (111.37
 mM)δ½ ε ζ»ζ §:ε °δε ι°He• η °ε -θ—δI ε• η η@Aδε-oO ι θΩ ε• 12ε• • ε• εI-εH'θ½»η ·ζ §ε• θ— θ ε ι 'η» θ• —oO ζ— θ
 °Y δI ε• ε ι • Ro 41-3696oO LirequiniloO oO 1γ • 3γ • 5ε 10 mgoO γ • zolpidemoO 10
 mgoO ε ε ζ °ε ε—HeΩ η• θΩ• ε °θ"η °ε θ °εΩ η ε½±ε • γ εΩ η• θΩ• ε °θ"η °ζ" θ— oO θ• θE ε ζ°ζ • ε δ½ δεIζ ε
 ε δ»»ε 'η °ζ±η δε ι °ε oO ε ε «ε °η °θ• —ε • γ • η °θ• —ε• 1.5ε° ζ Ae 8ε° ζ AθΩ θ' γ θ °εΩ ζ" θ— ε °η °θ• —ε• 8
 ε° ζ AθΩ θ' oO ζA ε• ε εΩ ε °η °θ• —ε• 1.5ε° ζ Ae-δH η 15δε ε• θ—• ε θ"γ η °θ• —ε• 1.5ε° ζ AoO 10 mg
 zolpidemεO θ"η εΩ η• θΩ• ε °ζ εI ζ °Yθ ε §δI Ro 41-3696η δ»δ½ ε ι • γ 5 mge 10 mg Ro
 41-3696δ» ε• zolpidemη ζ ζ ζ °Yθ ε §δI ε ζ °ε oO P<0.05oO γ ζ-ζ ζ @ζ °oO η °θ• —ε• 8ε° ζ AoO 5 mge 10 mg Ro
 41-3696η θ½»εY ζ η ζ εI δε ε ζ °ε η °Eζ— ζ η» θ 'ε-ζ • δH δε η ζ °Yθ ε ε «oO θ zolpidemη ζ ζ δε ε ζ
 °ε ζ ε• εO γ θ °εΩ ζ" θ— η» ζ ζ °Yη loO ε-δH ε• ε εΩ θ ½ε ε• ε• zolpidemζ ζ ζ °Yε½±ε • γ ε—HδI ι« θ°Y5
 mgn ε ι • oO ζ θ§ ε— ε °Ro 41-3696ε—HθΩ δI ε• ι • η ε½±ε • γ ζ »δH oO η»• ζ" θ— η ζ ζ ε ι • η Ro
 41-3696ε °η °θ• —ε• 1.5ε° ζ Ae—HeΩ η• θΩ• ε °θ"η °ε θ °εΩ η ε½±ε • ε• ε° δI 10 mg zolpidemoO δ½ Ro
 41-3696η ζ ζ δOOδH ζ • η»-ζ Ai 'ζ 'ι Ωγ [2]δΩ'ε• ι θ—:η!ε• ι ι • |η η»• η§ ε-ι• Aη °H:GABA Receptor

Note: This product is for in vitro research use only and is not intended for use in humans or animals.