

4-Methylherniarin

Catalog No: #S08814



Package Size: #S08814-1 1 mL #S08814-2 200 mg #S08814-3 500 mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	4-Methylherniarin
Other Names	Bacterial;7-Methoxy-4-methylcoumarin;4-Methylherniarin;4-Methylherniarin;4Methylherniarin;inhibit;Inhibitor;ηΥ η ²ι θ± η´ η ²ι
Calculated MW	190.2
Storage	Powder: -20°C for 3 years . In solvent: -80°C for 1 year . Shipping with blue ice.

Application Details

ε ε• εO• :C11H10O3Smile:COC1=CC=C2C(C)=CC(=O)OC2=C1DMLε• :MFCD00009773η¹—ε¹;97.31%ζ¹Aθ§εε¹;DMSO:11 mg/mL (57.83 mM)• Aη Ηζ´»ζ §:S. sonnei:13.47 ΞOg/ml (IC50)|B. subtilis:11.76 ΞOg/ml (IC50)δΩ'ε• ·ι θ—:εΥ η η ©ε-|• Aη Η:Antibacterial

Note: This product is for in vitro research use only and is not intended for use in humans or animals.