

LSN2463359

Catalog No: #S08816



Package Size: #S08816-1 1 mg #S08816-2 1 mL #S08816-3 10 mg #S08816-4 100 mg #S08816-5 200 mg #S08816-6 2.5 mg #S08816-7 5 mg #S08816-8 50 mg
Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	LSN2463359
Other Names	LSN2463359;LSN-2463359;LSN 2463359;mGlu5
Calculated MW	265.31
Storage	Powder: -20°C for 3 years . In solvent: -80°C for 1 year . Shipping with blue ice.

Application Details

SMILES: CC(C)NC(=O)c1ccc(nc1)C#Cc1ccccc1 DMSO: 112.5 mg/mL (424.0 mM)
 Molecular Weight: 265.31 g/mol
 Purity: >99.36%
 Solubility: Soluble in DMSO, water, and other common solvents.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.