

LX-6171

Catalog No: #S08819



Package Size: #S08819-1 1 mg #S08819-2 1 mL #S08819-3 10 mg #S08819-4 100 mg #S08819-5 200 mg #S08819-6 2.5 mg #S08819-7 5 mg #S08819-8 50 mg
Others: order@signalwayantibody.com Support: tech@signalwayantibody.com

Description

Product Name	LX-6171
Other Names	LX 6171;LX6171;LX-6171;SLC6A7
Calculated MW	377.87
Storage	Powder: -20°C for 3 years . In solvent: -80°C for 1 year . Shipping with blue ice.

Application Details

ε ε• εO• :C22H20CIN3OSmile:C1C=C(C2=CC=C(C(=O)C3CCN(CC3)C=4N=CC=CN4)C=C2)C=CC1η1—ε1;97.06%ζ|Aθ\$£ε1;DMSO:25 mg/mL (66.16 mM)δΩ'ε• ·1 θ·—:ε Aδ» ·· Aη 'H:Others

Note: This product is for in vitro research use only and is not intended for use in humans or animals.