

Recombinant Flag-Tag Monoclonal Antibody(3O2), HRP Conjugated

Catalog No: #T402

Package Size: #T402-1 50ul #T402-2 100ul



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Flag-Tag Monoclonal Antibody(3O2), HRP Conjugated
Purification	affinity purified by Protein A
Applications	WB,ELISA
Species Reactivity	Species independent
Conjugates	HRP
Formulation	1Γ PBS (pH7.4) oO 0.05% proclin300oO 50% Glycerol
Storage	Store at -20°C/1 year. Avoid repeated freeze/thaw cycles.

Application Details

WB 1:30000-1:100000

ELISA 1:5000-1:10000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.