

Recombinant His-Tag Monoclonal Antibody(14W7), HRP Conjugated

Catalog No: #T404

Package Size: #T404-1 50ul #T404-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant His-Tag Monoclonal Antibody(14W7), HRP Conjugated
Purification	affinity purified by Protein A
Applications	WB,ELISA
Species Reactivity	Species independent
Conjugates	HRP
Formulation	1Γ PBS (pH7.4) oO 0.05% proclin300oO 50% Glycerol
Storage	Store at -20°C/1 year. Avoid repeated freeze/thaw cycles.

Application Details

WB 1:30000-1:100000ELISA 1:5000-1:10000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.