

SRT-Tag Mouse Monoclonal Antibody(4H9)

Catalog No: #T518



Package Size: #T518-1 50ul #T518-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SRT-Tag Mouse Monoclonal Antibody(4H9)
Host Species	Mouse
Clonality	Monoclonal
Purification	Affinity purification using immunogen.
Applications	WB
Specificity	The SRT Tag Mouse Monoclonal Antibody detects SRT Tag recominant protein.
Target Name	SRT-Tag
Other Names	SRT Tag antibody
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 51% Glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:5000~1:10000

Images



Western blot analysis of SRT-Recombinant protein, using #T518 diluted at 1) 1:5,000 2) 1:10,000

Background

The SRT-tag is a peptide epitope comprised of 10 amino acids(TFIGAIATDT) derived from the crystalline surface layer protein of Rickettsia typhi1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.