

RFP-Tag Rabbit Polyclonal Antibody

Catalog No: #T522



Package Size: #T522-1 50ul #T522-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	RFP-Tag Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purification using immunogen.
Applications	WB
Species Reactivity	Human;Mouse;Rat
Specificity	The RFP tag antibody can recognize RFP tag fusion proteins.
Conjugates	Unconjugated
Target Name	RFP-Tag
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:5000~1:10000

Images



2ug RFP fusion protein+ Primary antibody #T522 dilution at 1:5,000 (lane 1) and 1:10,000 (lane 2).

Background

The red fluorescent protein cloned from *Discosoma* coral (DsRed or drFP583) holds great promise for biotechnology and cell biology as a spectrally distinct companion or substitute for the green fluorescent protein (GFP) from the *Aequorea* jellyfish.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.