

TAP-Tag Mouse Monoclonal Antibody

Catalog No: #T527

Package Size: #T527-1 50ul #T527-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TAP-Tag Mouse Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	Affinity purification using immunogen.
Applications	WB IP
Specificity	Antibody detects TAP recombinant protein.
Target Name	TAP-Tag
Concentration	1.0mg/ml
Formulation	Mouse IgG1 in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:5000~1:10000

Images



Western blot analysis of TAP recombinant protein, using #T527 diluted at 1:5,000.

Background

Epitope mapping within the calmodulin binding peptide (CBP) region of the TAP (tandem affinity purification) tag of human.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.