

RFP-Tag Mouse Monoclonal Antibody

Catalog No: #T608

Package Size: #T608 100ul



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	RFP-Tag Mouse Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Purification	Affinity-chromatography
Applications	WB IHC IF/ICC
Specificity	RFP-Tag mouse monoclonal antibody detects endogenous levels of internal, C-terminal, and N-terminal RFP recombinant proteins.
Immunogen Type	Recombinant protein
Immunogen Description	Full length RFP protein.
Target Name	RFP-Tag
SDS-PAGE MW	N/A
Concentration	1.0mg/ml
Formulation	Mouse IgG1 in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20 °C.Stable for 12 months from date of receipt


Application Details

WB: 1:3000~1:10000

IHC: 1:200

IF/ICC: 1:200

Images



1ug RFP fusion protein+ Primary antibody dilution at
1γ • 1:3,000 2γ • 1:5,000 3γ • 1:10,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.