

ECFP-Tag Monoclonal Antibody

Catalog No: #T616



Package Size: #T616-1 50ul #T616-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ECFP-Tag Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was affinity-purified from antiserum by affinity-chromatography
Applications	WB
Specificity	The antibody detects endogenous ECFP-Tag protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant Protein
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1:2000-5000

Images



Western blot analysis of Recombinant ECFP Protein using ECFP-Tag Monoclonal Antibody.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.