

Flag-Tag Rabbit mAb

Catalog No: #T617



Package Size: #T617-1 50ul #T617-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Flag-Tag Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB,ICC/IF,IP
Species Reactivity	Human,Mouse,Rat,All
Immunogen Description	A synthesized peptide
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.


Application Details

WB: 1:5000~1:50000

ICC/IF: 1:100~1:500

IP: 1:50~300

Images



All lanes use the Antibody at 1:2W dilution for 1 hour at room temperature.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.