

## Anti-HCoV-229E Spike Protein Recombinant Antibody

Catalog No: #V1002



Package Size: #V1002-1 50ug #V1002-2 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |                                                                                                                                                                           |
|--------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name       | Anti-HCoV-229E Spike Protein Recombinant Antibody                                                                                                                         |
| Host Species       | Human                                                                                                                                                                     |
| Clonality          | Monoclonal                                                                                                                                                                |
| Clone No.          | VS1322                                                                                                                                                                    |
| Isotype            | IgG1, kappa                                                                                                                                                               |
| Purification       | Protein A/G purified from cell culture supernatant.                                                                                                                       |
| Applications       | ELISA                                                                                                                                                                     |
| Species Reactivity | HCoV-229E                                                                                                                                                                 |
| Conjugates         | Unconjugated                                                                                                                                                              |
| Target Name        | E2, Peplomer protein, S, S glycoprotein, Spike glycoprotein                                                                                                               |
| Accession No.      | P15423                                                                                                                                                                    |
| Formulation        | 0.01M PBS pH 7.4                                                                                                                                                          |
| Storage            | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt. |

## Application Details

Purity:&gt;95% as determined by SDS-PAGE.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.