

Anti-PEDV Spike glycoprotein (RBD) Recombinant Antibody

Catalog No: #V1004

Package Size: #V1004-1 50ug #V1004-2 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Anti-PEDV Spike glycoprotein (RBD) Recombinant Antibody
Host Species	Human
Clonality	Monoclonal
Clone No.	VS1117
Isotype	VHH-Fc
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA
Species Reactivity	Porcine epidemic diarrhea virus (PEDV)
Conjugates	Unconjugated
Target Name	Spike glycoprotein, S glycoprotein, E2, Peplomer protein, S
Accession No.	Q91AV1
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

Purity:>95% as determined by SDS-PAGE.ELISA: 1:2000-1:20000

Images



SDS-PAGE for Anti-PEDV Spike glycoprotein (RBD) Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.