

Anti-SADS-CoV Spike Recombinant Antibody

Catalog No: #V1005



Package Size: #V1005-1 50ug #V1005-2 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Anti-SADS-CoV Spike Recombinant Antibody
Host Species	Human
Clonality	Monoclonal
Clone No.	VS2147
Isotype	IgG
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA, Neutralization
Species Reactivity	Porcine enteric alphacoronavirus GDS04
Conjugates	Unconjugated
Target Name	spike glycoprotein [Porcine enteric alphacoronavirus GDS04], S, N-terminal domain of the S1 subunit of coronavirus Spike (S) proteins; cl40465, SD-1 and SD-2 subdomains, the S1/S2 cleavage region, and the S2 fusion subunit of the CoV spike (S) glycoprotein from Rhinolophus bat coronavirus HKU2 and related alphacoronaviruses; cd22371, S1/S2 cleavage region, putative fusion peptide, internal fusion peptide, heptad repeat 1 [structural motif], heptad repeat 2 [structural motif]
Accession No.	ASK51717.1
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

Purity:>95% as determined by SDS-PAGE.ELISA: 1:2000-1:20000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.