

Anti-influenza A virus (H5N1) Non-structural protein 1
Recombinant Antibody

Catalog No: #V1006



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Package Size: #V1006-1 50ug #V1006-2 100ug


Description

Product Name	Anti-influenza A virus (H5N1) Non-structural protein 1 Recombinant Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	VS1341
Isotype	IgG1
Purification	Protein A/G purified from cell culture supernatant.
Applications	WB, ELISA
Species Reactivity	Influenza A virus (H5N1), Influenza A virus (H1N1)
Conjugates	Unconjugated
Target Name	NS1
Accession No.	O56264
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

Purity:>90% as determined by SDS-PAGE.ELISA: 1:2000-1:20000; WB: 1:500-1:2000

Images



SDS-PAGE for Anti-influenza A virus Non-structural protein 1 Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.