

Anti-Influenza B virus/Flu-B virus Monoclonal Antibody

Catalog No: #V1007




Package Size: #V1007-1 50ug #V1007-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Anti-Influenza B virus/Flu-B virus Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	VS1938
Isotype	IgG
Purification	Protein A/G purified from cell culture supernatant.
Applications	GICA
Species Reactivity	Influenza B virus (Flu-B)
Immunogen Description	Influenza B virus/Flu-B virus Protein.
Conjugates	Unconjugated
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Images



SDS-PAGE for Anti-Influenza B virus/Flu-B virus Monoclonal Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.