

Anti-Influenza A virus Neuraminidase/NA protein Recombinant Antibody

Catalog No: #V1010

Package Size: #V1010-1 50ug #V1010-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Anti-Influenza A virus Neuraminidase/NA protein Recombinant Antibody
Host Species	Human
Clonality	Monoclonal
Clone No.	VS2248
Isotype	IgG1
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA
Species Reactivity	Influenza A virus (strain A/Tern/Australia/G70C/1975 H11N9)
Conjugates	Unconjugated
Target Name	Neuraminidase, NA
Accession No.	P03472
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

Application Details

Purity:>95% as determined by SDS-PAGE.ELISA: 1:2000-1:20000

Images



SDS-PAGE for Anti-Influenza A virus Neuraminidase/NA protein Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.