

InVivoMAb Anti-Influenza A virus(strain A/Japan/305/1957 H2N2) Neuraminidase/NA Antibody

Catalog No: #V1011

Package Size: #V1011-1 100ug



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	InVivoMAb Anti-Influenza A virus(strain A/Japan/305/1957 H2N2) Neuraminidase/NA Antibody
Host Species	Human
Clonality	Monoclonal
Clone No.	VS2023
Isotype	IgG1
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA, Neutralization
Species Reactivity	Influenza A virus (strain A/Japan/305/1957 H2N2) (Q1K9Q1), Influenza A virus (strain A/RI/5-/1957 H2N2) (P03484)
Conjugates	Unconjugated
Target Name	Neuraminidase, 3.2.1.18, NA
Accession No.	Q1K9Q1
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

Purity:>95% as determined by SDS-PAGE.

Images



SDS-PAGE for InVivoMAb Anti-Influenza A virus Neuraminidase/NA Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.