

Anti-Influenza A virus Neuraminidase/NA Recombinant Antibody

Catalog No: #V1016

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #V1016-1 50ug #V1016-2 100ug


Description

Product Name	Anti-Influenza A virus Neuraminidase/NA Recombinant Antibody
Host Species	Human
Clonality	Monoclonal
Clone No.	VS1707
Isotype	IgG1
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA
Species Reactivity	Influenza A virus (A/environment/Shanghai/S1439/2013(H7N9)) (S5MF06), Influenza A virus (A/Shanghai/02/2013(H7N9)) (R4NFR6), Influenza A virus (A/chicken/Jiangsu/S1460/2016(H7N9)) (A0A2D0Z4C5)
Conjugates	Unconjugated
Target Name	Neuraminidase, 3.2.1.18, NA
Accession No.	S5MF06
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

Application Details

Purity:>95% as determined by SDS-PAGE.ELISA: 1:2000-1:20000

Images



SDS-PAGE for Anti-Influenza A virus Neuraminidase/NA Antibody .

Note: This product is for in vitro research use only and is not intended for use in humans or animals.