

Anti-Influenza A virus Neuraminidase/NA Recombinant Antibody

Catalog No: #V1021

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #V1021-1 50ug #V1021-2 100ug

Description

Product Name	Anti-Influenza A virus Neuraminidase/NA Recombinant Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	VS1537
Isotype	IgG1, kappa
Purification	Protein A/G purified from cell culture supernatant.
Applications	BIF, ELISA, Neutralization, X-ray crystallography
Species Reactivity	Influenza A virus (strain A/Japan/305/1957 H2N2)
Conjugates	Unconjugated
Accession No.	Q1K9Q1
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

Purity:>95% as determined by SDS-PAGE.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.