

Anti-Influenza A/B virus NA/Neuraminidase Broad-Neutralizing Recombinant Antibody

Catalog No: #V1022

Package Size: #V1022-1 50ug #V1022-2 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Anti-Influenza A/B virus NA/Neuraminidase Broad-Neutralizing Recombinant Antibody
Host Species	Human
Clonality	Monoclonal
Clone No.	VS1449
Isotype	IgG
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA, Neutralization
Species Reactivity	Influenza A virus, Influenza B virus
Conjugates	Unconjugated
Target Name	EC:3.2.1.18, NA, Neuraminidase
Accession No.	P03468
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

Purity:>95% as determined by SDS-PAGE.ELISA: 1:2000-1:20000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.