

Anti-Influenza A virus M/Matrix protein 1 Recombinant Antibody

Catalog No: #V1023

Package Size: #V1023-1 50ug #V1023-2 100ug



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Anti-Influenza A virus M/Matrix protein 1 Recombinant Antibody
Host Species	Human
Clonality	Monoclonal
Clone No.	VS1690
Isotype	IgG1
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA, BLI
Species Reactivity	Influenza A virus (H1N1)
Conjugates	Unconjugated
Target Name	M1
Accession No.	P03485
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

Purity:>90% as determined by SDS-PAGE.ELISA: 1:2000-1:20000

Images



SDS-PAGE for Anti-Influenza A virus M/Matrix protein 1 Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.