

Anti-Influenza A virus (H5N1) HA/Hemagglutinin Recombinant Antibody

Catalog No: #V1025



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Package Size: #V1025-1 50ug #V1025-2 100ug


Description

Product Name	Anti-Influenza A virus (H5N1) HA/Hemagglutinin Recombinant Antibody
Host Species	Human
Clonality	Monoclonal
Clone No.	VS2244
Isotype	IgG1, kappa
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA, IF, SPR, WB
Species Reactivity	Influenza A virus (A/Viet Nam/1203/2004(H5N1)), Influenza A virus (strain A/Goose/Guangdong/1/1996 H5N1 genotype Gs/Gd)
Conjugates	Unconjugated
Target Name	HA
Accession No.	Q6DQ33, Q9Q0U6
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

Application Details

Purity:>95% as determined by SDS-PAGE.ELISA: 1:2000-1:20000; WB: 1:500-1:2000

Images



SDS-PAGE for Anti-Influenza A virus (H5N1)
HA/Hemagglutinin Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.