

InVivoMAb Anti-Influenza B virus(strain B/Maryland/1959) HA/Hemagglutinin Antibody

Catalog No: #V1026

Package Size: #V1026-1 100ug



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	InVivoMAb Anti-Influenza B virus(strain B/Maryland/1959) HA/Hemagglutinin Antibody
Host Species	Human
Clonality	Monoclonal
Clone No.	VS1525
Isotype	IgG1, kappa
Purification	Protein A/G purified from cell culture supernatant.
Applications	Blocking, ELISA
Species Reactivity	Influenza B virus (strain B/Maryland/1959)
Conjugates	Unconjugated
Target Name	Hemagglutinin, Hemagglutinin HA1 chain, Hemagglutinin HA2 chain, HA
Accession No.	P03461
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

Purity:>95% as determined by SDS-PAGE.

Images



SDS-PAGE for InVivoMAb Anti-Influenza B virus HA/Hemagglutinin Antibody .

Note: This product is for in vitro research use only and is not intended for use in humans or animals.