

InVivoMAb Anti-Influenza A
virus(A/Perth/16/2009(H3N2)) HA/Hemagglutinin
Antibody

Catalog No: #V1027

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Package Size: #V1027-1 100ug


Description

Product Name	InVivoMAb Anti-Influenza A virus(A/Perth/16/2009(H3N2)) HA/Hemagglutinin Antibody
Host Species	Human
Clonality	Monoclonal
Clone No.	VS1055
Isotype	IgG1, kappa
Purification	Protein A/G purified from cell culture supernatant.
Applications	Blocking, ELISA
Species Reactivity	Influenza A virus (A/Perth/16/2009(H3N2))
Conjugates	Unconjugated
Target Name	Hemagglutinin, Hemagglutinin HA1 chain, Hemagglutinin HA2 chain, HA
Accession No.	C6KNH7
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

Purity:>95% as determined by SDS-PAGE.

Images



SDS-PAGE for InVivoMAb Anti-Influenza A virus
HA/Hemagglutinin Antibody .

Note: This product is for in vitro research use only and is not intended for use in humans or animals.