

Anti-Influenza A Virus H7 subtype Hemagglutinin/HA Recombinant Antibody

Catalog No: #V1030



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Package Size: #V1030-1 50ug #V1030-2 100ug


Description

Product Name	Anti-Influenza A Virus H7 subtype Hemagglutinin/HA Recombinant Antibody
Host Species	Human
Clonality	Monoclonal
Clone No.	VS1185
Isotype	IgG1, kappa
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA, IF, BLI, in vivo protection, in vitro neutralization
Species Reactivity	Influenza A virus (H7N1), Influenza A virus (H7N3)
Conjugates	Unconjugated
Target Name	Hemagglutinin (H7 subtype)
Accession No.	P03459, P03458
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

Purity:>90% as determined by SDS-PAGE.ELISA: 1:2000-1:20000

Images



SDS-PAGE for Anti-Influenza A Virus H7 subtype Hemagglutinin/HA Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.